
Kaltlichtquelle Coldlightsource

intralux[®] dc - 1100

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

Unit and its components

A. Starting up

1. First check, whether the factory - set mains corresponds with the local mains voltage.
2. Loosen cover screws (9) on both sides and remove cover (8).
3. Connect halogen lamp (4) into socket (3) and push it down into the lamp holder to a positive - stop. Do not touch the bulb or the reflecting surface of the lamp!
4. Replace cover (8) and refasten cover screws (9).
5. Plug the mains supply cable into the power entry unit (10) and connect to the mains.
6. Switch throw - over switch (14) at the rear housing in position „Local“. (Only for units with the serial interface!)
7. Switch the power switch (7) on the front of the instrument to position 1.

Caution:

- Installation and maintenance only by qualified persons.
- The safety of the lightsource can only be guaranteed if used according to the operation instructions.
- Clean the lightsource with a moist cloth. Unplug the lightsource first !
- Lamp life depends to a large extent on its operating temperature. To achieve best results, please observe that the air can flow freely (Air inlet at bottom of housing, air outlet on the rear of housing)!

B. Operation

1. Plug the mains supply cable into the power entry unit (10) and connect to the mains.
2. Switch throw - over switch (14) at the rear housing in position „Local“. (Only for units with the serial interface!)
3. Switch the power switch (7) on the front of the instrument to position 1.
4. Adjust intensity with the intensity control knob (6). (This regulation of intensity changes the colour temperature of the lamp and affects his life.)
5. When the intensity setting has been selected the light supply can be adjusted continuously by means of the brightness control (5).

Notice:

Between the jacks (for banana plugs) at the rear housing - below the fan - the actual lamp voltage can be measured.

C. Lamp replacement

Caution:

- Lamp should be cooled off before touching!

1. Set the power switch (7) on the front of the instrument to position **0** and disconnect the mains supply cable from the instrument and from the mains.
2. Remove cover according to A.2.
3. Remove the defective lamp (4) together with the lamp socket (3) from the lamp holder and remove it from the lamp socket (3). Replace the lamp (4) with a new one (EKE or DDL - type).
4. Further steps as per A.3. – A.7.

D. Fuse replacement

1. Set the power switch (7) on the front of the instrument to position **0** and disconnect the mains supply cable from the instrument and from the mains.
2. Pull out fuse holder (11) on the rear of the housing.
3. Replace the defective fuses (2) with spare fuses. (Replace only of same type as per manufacturer's label!)
4. Close the fuse holder (11).
5. Further steps as per A.5. – A.7.

E. Lamp socket replacement

Caution:

- If the lamp socket (3) shows signs of burning and wear between the lamp pin contacts it must be replaced.
- Lamp should be cooled off before touching!

1. Set the power switch (7) on the front of the instrument to position **0** and disconnect the mains supply cable from the instrument and from the mains.
2. Remove cover according to A.2.
3. Remove the halogen lamp (4) together with the lamp socket (3) from the lamp holder and remove it from the lamp socket (3). Do not touch the bulb or the reflecting surface of the lamp!
4. Separate the old or damaged lamp socket (3) from the cable terminal (12) and replace it with a new one.
5. Further steps as per A.3. – A.7.

F. Serial interface RS-232

(optional)

Please find the description of the serial interface in the rear part of this manual.

Spare Parts

1. Heat protection filter	91611.006
2. Fuse (2.5A T / 250V)	90045.118
3. Lamp socket.....	25750.045
4. Halogen lamp EKE 21V / 150W	90357.006
Halogen lamp DDL 20V / 150W	23500.014
– Mains supply cable Europe	90032.111
Mains supply cable Switzerland	90032.112
Mains supply cable GB.....	90032.117

Legend

5. Brightness control knob	(crescent shaped diaphragm)
6. Intensity control knob	
7. Power switch	
8. Cover	
9. Cover screw	
10. Power entry unit	
11. Fuse holder	
12. Cable terminal	
13. Serial interface RS-232	(optional)
14. Throw - over switch	(optional)

Ambient conditions

For operating:

Ambient temperature	0°C to 40°C
Relative humidity (RH) (non condensing)	5% to 95%
Overvoltage category	class II
Pollution degree	2

For storage:

Ambient temperature	-25°C to +70°C
Relative humidity (RH)	5% to 95%


Electrical data

Line Voltage	230VAC / 50Hz / 60Hz
Maximum power consumption	230W
Fuse	2.5AT / 250V (slow blow)

Serial Interface RS-232 (optional)

Cable:

The connection of the computer and the dc-1100s lightsource need two wires: Signal and Ground.

Computer / DB9F (frontview of plug) Lightsource / DB9M (frontview of plug)

Setup of the serial interface:

The serial interface needs following setup:

Baud: 9600	Databits: 7	Stopbit: 1	Parity: even
------------	-------------	------------	--------------

Communication:

The interface contains a 10-bit DAC. That means that the range of brightness can be adjusted from 0 (dark) to 1023 (max. intensity). The DAC features a serial input. To program the desired brightness 12 bits must be transmitted. First, 10 bits data (MSB first), second, 2 dummy bits. For each bit a whole RS-232 data package has to be transmitted. To send a 0, the value of this package is 0 too, to send a 1 the value is 255 (0xff).

Due to the internal structure of the interface the value of the brightness has to be „inverted“, i.e. subtracted from 1023.

Example: Program the brightness value to 555

1. Step: invert the value
 $1023 - 555 = 468$
2. Step: convert 468 in a 10-digit binary number
 $468 \text{ (dec.)} = 0'1111'010'100 \text{ (bin.)}$
3. Step: transmit 10 bits
send(0), send(255), send(255), send(255), send(0), send(255),
send(0), send(255), send(0), send(0).
4. Step: transmit 2 dummy bits
Send(0), send(0).

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Referee also to the procedure „serout“ in the program „serial.c“ on the enclosed floppydisk.

```
void serout(short inten, short prt)
{
    short j=512,i;
    short adr=COM1;                /* use com1 */

    if (prt==2) adr=COM2;          /* use com2 */
    for (i=0;i<10;i++)            /* decut the value of the intensity */
    {                               /* to its binary components and */
        if (inten&j) outportb(adr,0xff); /* send it bit by bit... */
        else outportb(adr,0x00);        /* send 0xff for 1, 0x00 for 0 */
        delay(2); /* wait unit done */ /* send MSB first */
        j/=2;
    }
    outportb(adr,0x00);            /* send 2 dummy bits due to */
    delay(2);                      /* internal structure of the DAC */
    outportb(adr,0x00);
    delay(40);                      /* wait 40ms to recover */
}
```